SMARTBRANE
Resorbable Pericardium Membrane

SMARTBRANE

Resorbable Pericardium Membrane

LV	4	Lietošanas instrukcija SMARTBRANE
ET	9	SMARTBRANE'i kasutusjuhend
DE	14	Gebrauchsanweisung für SMARTBRANE
EN	19	Instructions for use for SMARTBRANE
FR	24	Mode d'emploi pour SMARTBRANE
IT	29	Istruzioni per l'uso di SMARTBRANE
ES	34	Instrucciones de uso para SMARTBRANE
SV	39	Bruksanvisning för SMARTBRANE
NL	44	Instructies voor het gebruik van SMARTBRANE

Latviešu

Produkta apraksts

Materiāls

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Product Description

Material

SMARTBRANE is a natural implantable resorbable collagen membrane for use in periodontal and/or dental surgery procedures to support guided tissue regeneration (GTR) and guided bone regeneration (GBR). The membrane is manufactured using a controlled multistage purification process. SMARTBRANE is manufactured from porcine pericardium from veterinary controlled pigs from the Netherlands. The pericardium are thoroughly purified, degreased, lyophilized and sterilized in a double peelable pouch using gamma sterilization.

Properties

SMARTBRANE functions as a barrier between bone graft material and soft tissue. The membrane serves as a resorbable scaffold that is eventually resorbed and replaced by tissue.

Product forms

SMARTBRANE is available in different sizes.

Intended use

SMARTBRANE is intended to aid in the regeneration and integration of oral tissue components in guided bone regeneration (GBR) and guided tissue regeneration (GTR) procedures.

Indications

SMARTBRANE is intended for use in intraoral/maxillofacial surgery procedures as a material for placement in the area of periodontal defects, dental implants, bone defects or ridge reconstructions to aid in wound healing post-surgery.

Contraindications

- Known hypersensitivity to (porcine) collagen.
- Active or latent infection at or around the wound site.
The practitioner should be confident that any active or recent infection has been properly treated.
- Any systemic disorder or disease that involves an unacceptable increase in postoperative risks or complications.

Warnings

Pregnant women

SMARTBRANE has not been tested in pregnant women.

Nursing mothers

It is not known whether degradation products of SMARTBRANE are secreted in human milk.

Pediatric use

SMARTBRANE is NOT intended for use in pediatric patients.

Caution/Precautions

- SMARTBRANE is designed for single patient use. Do NOT

re-use or re-sterilize. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury or serious illness of the patient. Unused portions must be discarded and not be reused.

- SMARTBRANE must NOT be used if the sterile package is opened or damaged.
- Do not use the membrane unless package integrity can be verified, proper storage has been assured and the expiration date has been verified. If the product, packaging or labeling is compromised in any way, do not use the membrane and notify your supplier.
- Store the sterile membrane according to the instructions on the label in a clean, dry place protected from sunlight.
- When handling the membrane, use sterile talc-free gloves to avoid contamination with particulates.
- Improper storage or handling may damage and/or contaminate the membrane and affect the sterility, handling characteristics, resorption characteristics and/or physical integrity of the device.
- As with any surgical procedure, infection is a risk. Use sterile intra-operative technique.
- The following complications are common with the surgical intervention with barrier membranes and therefore may not be totally excluded: Soft tissue dehiscence, hematoma, pain, inflammation, increased sensitivity and redness.

Procedure

Pre-operative

- Screen patients for a history of hypersensitivity for (porcine) collagen prior to implantation.
- Do not remove the membrane from the inner sterile pouch until just prior to implantation. Avoid contact between the membrane and surgical drapes and/or cotton gauze to avoid contamination with fibers.

Intra-operative

- A strict adherence to sterile protocol has to be followed du-

ring the implantation of the membrane.

- If the membrane is contaminated during the surgical procedure, it must be discarded.
- When the bone defect is exposed, it can be filled with a suitable material (e.g. allograft, autogenous, xenogenic or synthetic bone) if indicated.
- The membrane can be cut to size with a pair of sterile scissors to overlap the defect walls by at least 2 – 3 mm.
- After sizing, the membrane can, if preferred, be hydrated in a 0.9% sterile saline solution until it becomes soft and flexible.
- The membrane must be pressed down slightly after placement over the defect to ensure proper adhesion.
- The membrane can be sutured or fixed with nails or pins to prevent displacement.
- The mucoperiosteal flap must be repositioned over the membrane and sutured, completely covering the membrane.

Post-operative

- Patients should be monitored closely and instructed to inform the surgeon in case of adverse reactions.
- Patients should be informed to refrain from any activity that could cause the disturbance of the membrane.
- Patients should be instructed on a proper oral hygiene after surgery.

Treatment risks

- SMARTBRANE is solely an aid to bone grafting. Use of the device does not ensure successful treatment. The outcome of treatment depends on operative technique and host response.
- Infection is a risk. Use sterile intraoperative technique. Avoid use at infected sites.
- SMARTBRANE is NOT self-supporting and must be used in combination with a bone graft material in order to maintain space under the membrane.
- In rare cases, a local reaction can occur in patients without a

history of hypersensitivity to collagen. The following side effects may be observed: dehiscence with early membrane exposure, hematoma, swelling, increased sensitivity or pain, local redness and/or inflammation.

Storage and stability

SMARTBRANE must be kept in its original packaging, stored at room temperature in a dry place. Long term storage is recommended at room temperature in a dry place and protected from direct sunlight. Exposure to temperatures below -20°C or above 40°C for periods longer than 48 hours should be avoided as product stability has not been tested outside these temperature conditions. The expiration date is printed on the outer label. Do NOT use after expiration date.

Please note

- Practitioners must be trained in surgical techniques for oral tissue regeneration using a barrier membrane. Specialized training is strongly recommended.
- Practitioners must have appropriate knowledge and instruction in the handling of SMARTBRANE for using the product safely and properly in accordance with these instructions for use.
- SMARTBRANE must be used in accordance with the instructions for use provided by the manufacturer. It is the practitioner's responsibility to use the device in accordance with these instructions for use and to determine, if the device fits to the individual patient situation.

Validity

Upon publication of these instructions for use, all previous versions are superseded. ©REGEDENT AG. All rights reserved.

Availability

SMARTBRANE is not available in all countries.

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